

professional, glass body

PLANE pH ELECTRODE

Model : PE-08

ISO-9001, CE, IEC1010



↑
Plane head

FEATURES

- * Professional plane pH electrode, application for Laboratory, high quality & field usage, such as s skin, food... pH measurement.
- * 0 to 14 pH.
- * Glass body with BNC connector
- * Size : 15 mm Dia. x 137.5 mm.



Lutron

LUTRON ELECTRONIC

The Art of Measurement

PLATE

pH ELECTRODE

Model : PE-08

Features	General purpose, laboratory & field usage.
Measuring Range	1 to 13 pH (typical 0 to 14 pH).
Operation Temperature	5 to 80 °C (41 to 176 °F).
Electrode Structure	Combination type.
Zero Potential for pH Value	7 ± 1 pH.
Body Material	Glass
Connector	BNC.
Mechanical Protection	With protection bottle on the electrode head.
Dimensions	Body length - 140 mm. Body Dia. - 150 mm. Cable length - 1000 mm.
Applications	Professional plane pH electrode, application for Laboratory, high quality & field usage, such as s skin, paper, fruit, food, meat... pH measurement.

1408-PE08