

PH ELECTRODE (GENERAL PURPOSE) MODEL: EH-801

FEATURES

High quality, professional, laboratory & field usage.

SPECIFICATIONS

Measuring Range	0.00 to 14.00 pH	
Operating Temperature	5 to 65°C (41 to 149°F)	
Electrode Structure	Combination Type	
Reference of pH Electrode	Ag / AgCl, Double Junction	
Body Material	Epoxy	
Connector	BNC	
pH Value of Zero mV	7.00 ± 1 pH	
Repeatability	≥ 0.05 pH	
Mechanical Protection	With protection bottle on the electrode head	
Weight	68 grams (including cable & socket)	
Dimensions	Body Length	120 mm
	Body Diameter	10 mm
	Cable Length	1200 mm

