PEN pH METER

Model: PH-223 *ISO-9001, CE, IEC1010*







The Art of Measurement

build ATC sensor, economic type

PEN PH METER

Model: PH-223

FEATURES

* Pen type digital pH meter, all in one, pH electrode is		
included, easy for general purpose application.		
* Build in temperature sensor for auto temperature		
compensation (ATC).		
` pH 4, pH 7 standard solution are standard accessories.		
* Easy to change the pH electrode.		
* IP67, water proof and protection.		
* Auto calibration for pH 4, pH 7 or pH 10.		
* High input impedance.		
* LCD with two displays.		
* Data hold function for freezing the desired value.		
* Auto power off.		
* Records max. and min. value with recall.		
* Microcomputer circuit, intelligent function, high accuracy.		
* Compact size, light weight.		
* Power supply by DC 1.5 V battery (UM4/AAA) x 4 PCs,		
* Available for wide applications, such as aquarium,		
beverage, fish hatcheries, food processing,		
photography, laboratory, quality control, school &		
colleges, swimming pools, water conditions.		

SPECIFICATIONS

Display	LCD, size: 20 mm x 28 mm.		
	* With b	ar graph indicator.	
Measurement	рН	0 to 14 pH.	
Range			
Resolution	рН	0.01 pH.	
Accuracy	рН	± 0.02 pH	
		* After calibration, meter	
		without PE-12 electrode.	
		± 0.2 pH	
		* After calibration, meter	
		with PE-12 electrode.	
pH Electrode	Included. * PE-12		
pH Temp.	build thermister Temp. sensor for		
Compensation	automatic temperature compensation		
	(ATC).		

pH Input	10^12 ohms.	
Impedan:		
Data Hold	Freeze the display reading.	
Memory Recall	Maximum & Minimum value.	
Sampling	Approx. 0.8 second.	
Time		
рН	pH 7, pH 4 or pH 10, 3 points	
Calibration	calibration ensure the best linearity and	
	accuracy.	
Circuit	Custom one-chip of microprocessor LSI	
	circuit.	
Operating	Meter:	
Temperature	0 to 60 $^{\circ}\mathrm{C}$ (32 to 140 $^{\circ}\mathrm{F}$).	
	pH electrode :	
	0 to 80 $^{\circ}\mathrm{C}$ (23 to 176 $^{\circ}\mathrm{F}$).	
Operating	Less than 80% RH.	
Humidity		
Power Supply	DC 1.5V battery (UM-4/AAA) x 4 PCs.	
Power	Approx. 4.8 mA.	
Consumption		
Dimension	186 x 40 x 32 mm	
	(7.3 x 1.6 x 1.3 inch).	
	* Meter with pH electrode.	
Weight	131 g/0.36 LB (included electrode).	
Power off	Auto power shut off to save battery life.	
Standard	Instruction Manual 1 PC	
Accessories	pH electrode, PE-12 1 PC	
	pH 7.00 buffer solution, PH-07A 1 PC	
	pH 4.00 buffer solution, PH-04A 1 PC	
	* PE-12, PH-07A, PH-04A are the	
	consumer accessory.	
Optional	pH ELECTRODE, Model : PE-12	
Accessories	pH 7.00 buffer solution, Model : PH-07A	
	pH 4.00 buffer solution, Model : PH-04A	
	Hard carrying case, Model : CA-06	
without notice	Soft carrying case, Model : CA-05A	

1303-PH223

^{*} Appearance and specifications listed in this brochure are subject to change without notice.